

 <p>एफएसएसएआई fssai भारतीय खाद्य सुरक्षा और मानक प्राधिकरण Food Safety and Standards Authority of India स्वास्थ्य और परिवार कल्याण मंत्रालय Ministry of Health and Family Welfare</p>	Method for Determination of Iron in Vitamin Mineral Premix for Preparation of Fortified Rice Kernel (FRK)		
Method No.	FSSAI.VMP-FRK.16.008.2023	Revision No. & Date	0.0
Scope	<p>The Scope of this Method is Applicable for Quantification of Iron in Premix at 5000 mg/Kg LOQ Level (with respect to the Sample) by Using Atomic Absorption Spectroscopy (AAS).</p> <p>Limit of Detection 12.5 mg/Kg with respect to the Standard. Limit of Quantification 25.0 mg/Kg with respect to the Standard. Limit of Quantification 5000 mg/Kg with respect to the Sample.</p>		
Caution (Safety & Precautions)	<p>1. Concentrated Nitric Acid It is a Chemical which is corrosive to Metals. It causes severe skin burns and eye damage. It is toxic if inhaled. It is corrosive to the respiratory tract</p> <p>Following safety measures need to be taken during Handling of Concentrated Nitric Acid:</p> <ol style="list-style-type: none"> a) Do not breathe dust/fume/gas/mist/vapors/spray b) Wash face, hands and any exposed skin thoroughly after handling c) Wear protective gloves/protective clothing/eye protection/face protection d) Use only outdoors or in a well-ventilated area Keep away from heat/sparks/open flames/hot surfaces. e) Keep/Store away from clothing/ other combustible materials f) Take any precaution to avoid mixing with combustibles g) Keep only in original container h) Wear respiratory protection <p>2. Hydrogen Peroxide It is Oxidizing, Corrosive and Irritant chemical. Safety measures need to be taken during Handling of Hydrogen Peroxide: When handling moderate-to-high concentrations of Hydrogen Peroxide in the workplace, ensure eyewash stations and safety showers are accessible, and use splash goggles, gloves, and an approved Vapor Respirator.</p>		
Principle	Nitric acid, and hydrogen peroxide are added to homogenized Vitamin Premix sample in microwave vessels, and digested using a preprogramed temperature control. Analysis is performed by AAS.		
Apparatus/Instruments	<ol style="list-style-type: none"> 1. Atomic Absorption Spectroscopy - AAS 2. Microwave Digester 3. Analytical Balance 4. Micro Pipettes (20 -200 µL) & (100 -1000 µL) 		
Materials and Reagents	<ol style="list-style-type: none"> 1. Concentrated Nitric Acid (Purity- 69%) - Suprapure 2. Hydrogen Peroxide (Purity -30%) – LR Grade 3. CRM / Standard Stock Solution - Iron (Purity - 1000 mg/Kg) 		

	4. Purity of Argon and other gas, if used must fulfill the standard of instrument requirement																																										
Sample Preparation	<p><u>PREPARATION OF SAMPLE SOLUTION</u></p> <ol style="list-style-type: none"> 1. Weigh 0.50 g (\pm 0.05 g) of Homogenized Sample. 2. Transfer to Microwave Digestion Closed (MDC) Vessel. 3. Heated Milli Q Water at 60 °C. 4. Add 2.0 mL of Hot Milli-Q water. 5. Add 1.0 mL Hydrogen Peroxide. 6. Add 5.0 mL of Nitric Acid. 7. Close the Microwave Vessel tightly. 8. Kept at Room Temperature for 5 minutes. 9. Kept the Vessel rotor in Microwave Digester. 10. Cool the Vessel at Room Temperature after Digestion. 11. Add 10 mL of Milli Q water. 12. Mixed well. 13. Transfer to 100 mL Volumetric Flask. 14. Volume make-up to 100 mL with Milli-Q water. 15. Filter and use this for injecting on AAS. <p>Note: If required, dilute the sample for the desired concentration.</p>																																										
Method of Analysis (a) Preparation of Standard solutions	<p>A) <u>PREPARATION OF BLANK (5% NITRIC ACID)</u></p> <p>Transfer 7.25 mL of Nitric Acid (69%) in 100 mL Milli Q Water in Glass Bottle Mix well. Shake Vigorously.</p> <p>B) <u>PREPARATION OF CALIBRATION STANDARD SOLUTIONS</u></p> <p>Use Intermediate Standard Solution-1 for Preparing Calibration Standard Solutions as mentioned in below Table.</p> <table border="1"> <thead> <tr> <th>CAL. STANDARD SOLUTIONS</th> <th>SSS (mg/Kg)</th> <th>VOL. OF SSS (mL)</th> <th>VOL. OF NITRIC ACID (mL)</th> <th>FINAL VOL. (mL)</th> <th>FINAL CONC. (mg/Kg)</th> </tr> </thead> <tbody> <tr> <td>LS 6</td> <td>1000</td> <td>1.50</td> <td>0.5</td> <td>10</td> <td>150</td> </tr> <tr> <td>LS 5</td> <td>1000</td> <td>1.25</td> <td>0.5</td> <td>10</td> <td>125</td> </tr> <tr> <td>LS 4</td> <td>1000</td> <td>1.00</td> <td>0.5</td> <td>10</td> <td>100</td> </tr> <tr> <td>LS 3</td> <td>1000</td> <td>0.75</td> <td>0.5</td> <td>10</td> <td>75</td> </tr> <tr> <td>LS 2</td> <td>1000</td> <td>0.50</td> <td>0.5</td> <td>10</td> <td>50</td> </tr> <tr> <td>LS 1</td> <td>1000</td> <td>0.25</td> <td>0.5</td> <td>10</td> <td>25</td> </tr> </tbody> </table> <p>CAL : Calibration SSS : Standard Stock Solution VOL : Volume LS : Linearity Solution</p> <p>NOTE: Use Freshly Prepared Standard solutions for the Analysis.</p> <p>C) <u>PREPARATION OF BRACKETING STANDARD SOLUTION (50 mg/Kg)</u></p> <ol style="list-style-type: none"> 1. Transfer 0.5 ml from Standard Stock Solution of Iron (1000 mg/L) in 10 ml volumetric flask. 2. Add 0.5 ml nitric acid and made up the Volume till 10 ml volumetric flask by Milli-Q water and mix by Vortex Shaker Mixer. 	CAL. STANDARD SOLUTIONS	SSS (mg/Kg)	VOL. OF SSS (mL)	VOL. OF NITRIC ACID (mL)	FINAL VOL. (mL)	FINAL CONC. (mg/Kg)	LS 6	1000	1.50	0.5	10	150	LS 5	1000	1.25	0.5	10	125	LS 4	1000	1.00	0.5	10	100	LS 3	1000	0.75	0.5	10	75	LS 2	1000	0.50	0.5	10	50	LS 1	1000	0.25	0.5	10	25
CAL. STANDARD SOLUTIONS	SSS (mg/Kg)	VOL. OF SSS (mL)	VOL. OF NITRIC ACID (mL)	FINAL VOL. (mL)	FINAL CONC. (mg/Kg)																																						
LS 6	1000	1.50	0.5	10	150																																						
LS 5	1000	1.25	0.5	10	125																																						
LS 4	1000	1.00	0.5	10	100																																						
LS 3	1000	0.75	0.5	10	75																																						
LS 2	1000	0.50	0.5	10	50																																						
LS 1	1000	0.25	0.5	10	25																																						

(b) Instrument Details	a) Instrument : Atomic Absorption Spectrometer			
	b) Equipment Conditions : As detailed in below Table			
	Hollow cathode Lamp	Iron (Fe)		
	Lamp Current (mA)	5.0		
	Absorption Wavelength (nm)	372.0		
	Slit Width (nm)	0.2		
	Signal-Type	Atomic Absorption		
	Signal -Measurement	Integration		
	Oxidant	Air		
	Oxidant Flow (L/min)	13.5		
	Acetylene Flow (L/min)	2		
	Equation	Linear		
	Read Parameter			
	Time (sec)	10		
	Delay time (sec)	10		
c) Microwave Digestion Program				
SL. NO	RAMPING STAGE	HOLD TIME (Minutes)	TEMP (°C)	POWER (Watt)
1	1	20	180	800
2	2	10	160	800
3	3	10	140	800
4	COOL DOWN	10	-	-
Note: The make & model of Instrument can be changed. However, the Instrument should be able to achieve the desired LOD & LOQ value.				

Batch Organization	<u>Injection Sequence</u>		
	SL.NO.	NAME OF INJECTIONS	NUMBER OF INJECTIONS
	1	Blank	2
	2	Linearity Solution (LS) - 1	1
	3	Linearity Solution (LS) - 2	1
	4	Linearity Solution (LS) - 3	1
	5	Linearity Solution (LS) - 4	1
	6	Linearity Solution (LS) - 5	1
	7	Linearity Solution (LS) - 6	1
	8	Blank	2
	9	Sample Solution	1
	10	Blank	2
11	Bracketing Standard Solution	1	
TOTAL INJECTIONS		14	
Calculation with Units of Expression	<p>Carry out analysis and calculate Regression coefficient (R²) by analyzing the calibration standards by fitting the data into a linear regression curve, including zero.</p> <p>Calculate the Iron Content in Vitamin Premix using the following equation:</p> $\text{Iron (mg/Kg)} = \frac{\text{Instrument Conc. (mg/Kg)} \times \text{Make-up Volume (mL)}}{\text{Sample Weight (gm)}}$		
LOD & LOQ	<p>Limit of Detection 12.5 mg/Kg with respect to the Standard. Limit of Quantification 25.0 mg/Kg with respect to the Standard. Limit of Quantification 5000 mg/Kg with respect to the Sample.</p>		
Inference (Qualitative Analysis)	<p>PRT/MT/PRM/2023/001, Method Validation Protocol for Estimation of Iron in Premix Using Atomic Absorption Spectroscopy.</p> <p>AOAC 2011.14: Determination of Minerals and Trace elements in Milk & Milk Products, Infant Formula, and Adult Nutrition.</p>		
Approved by	Scientific Panel on Methods of Sampling and Analysis		